

Which Royal Blue tube should I use?

BD vacutainer, trace element, EDTA, plastic tube with lavender vertical strip on label

This tube is for:

Arsenic, Cadmium, Heavy Metals, Mercury, Cobalt WHOLE BLOOD, Chromium WHOLE BLOOD



Special order tube (order through Courier Supply Requisition)

BD vacutainer, trace element, no additive, serum plastic tube with red vertical strip on label

This tube is for:

Copper, Selenium, Iodine, Manganese, Zinc, Cobalt SERUM, Chromium SERUM



Available through SH Supply Chain ITM-1066808

Trace Element, No Additive, Aluminum Only

This tube is for:
Aluminum

Manufactures change frequently.

Please check tube label and lab catalog for correct tube. If you are still unsure, contact the Referral's Department.



Special order tube (order through Courier Supply Requisition)

IF YOU ARE UNSURE, CONTACT THE DEPARTMENT WHERE THE TEST IS PERFORMED