

SUPPLY REQUISITION

| SS | Orde | red by | by Phone () | |
|--|----------------------|---------------|---|-------|
| ORIGINAL WITH YOUR COURIER OR FAX TO 616.267. | | | | |
| PHLEBOTOMY | | QTY | REQUISITIONS, FORMS AND PATIENT EDUCA | |
| Tube gold 5 mL (SST) | (order as ea) | | Laboratory supply requisition | (orde |
| Tube red 6 mL (plain) | (order as ea) | | General lab/Custom requisition form | (10 |
| Tube lavender 3 mL (EDTA) | (order as ea) | | Anatomical pathology requisition form | (10 |
| Tube light blue 2.7 mL (3.2% citrate) | (order as ea) | | Advanced technology (ATL) requisition form | (10 |
| Tube dark green 6 mL (sod heparin) | (order as ea) | | Infectious disease (Microbiology) requisition form | (10 |
| Tube light green 4.5 mL (lith heparin PST) | (order as ea) | | Allergen test requisition form | (10 |
| Tube royal blue 6 mL (trace element - Red stripe) | (order as ea) | | Surgical Pathology requisition | (1 |
| Tube royal blue 6 mL MAYO (metal free EDTA - purple stripe) | (order as ea) | | Additional test requisition (add on) | (orde |
| Tube royal blue 6 mL MAYO (no additive aluminum) (limit 1) | (order as ea) | | Patient home collection instructions (English) | (orde |
| Tube tan 3 mL (EDTA) | (order as ea) | | Patient home collection instructions (Spanish) | (orde |
| Tube pink 6 mL (EDTA) | (order as ea) | | Orange Epic/Epic Care link lab order slips | (orde |
| Tube gray 6 mL (sod fluoride) | (order as ea) | QTY | CYTOLOGY | |
| Tube yellow 6 mL ACD (sol A) | (order as ea) | | Pap test thin prep vial w/ cervical broom | (|
| Tube yellow 6 mL ACD (sol B) | (order as ea) | | Pap test thin prep vial w/ cervical brush & spatula | (|
| Tube microtainer green (lith heparin) | (25/bag) | | Cytolyt for non GYN small mouth | (orde |
| Tube microtainer gold (SST) | (25/bag) | | Cytolyt for non GYN large mouth | (orde |
| Tube microtainer gold (SST) | (25/bag) (25/bag) | | Green top slide fixative | (orde |
| Tube microtainer levender (EDTA) | (25/bag) (25/bag) | | Fine needle aspirate (FNA) kit | , |
| · · | (| | | (orde |
| Tube QuantiFeron | (order as ea) | | UroVysion kits | (orde |
| Tube urine myoglobin (mayo) | (order as ea) | QTY | GLUCOLA | |
| Multi draw vacutainer needle 21 - 1¼" | (48/box) | | Glucola orange (50 gram) | |
| Multi draw vacutainer needle 22 - 1¼" | (48/box) | | Glucola orange (100 gram) | |
| Vacutainer holder | (each) | | Glucola fruit punch (100 gram) | |
| Clear aliquot tube | (25/bag) | | Glucola lemon lime (50 gram) | |
| Brown/Amber aliquot tube | (25/bag) | QTY | SPECIMEN BAGS | |
| MAYO metal free aliquot tube | (25/bag) | | White specimen bag (ambient) | (10 |
| MAYO clear aliquot tube | (100/bag) | | Green specimen bag (refrigerated) | (10 |
| TOXICOLOGY | | | Blue specimen bag (frozen) | (10 |
| Blue top drug screen cups | (order as ea) | | Red specimen bag (stat) | (1 |
| Drug screen cup seals (numbered) | (order as ea) | | Purple specimen bag (priority) | (10 |
| Drug screen seals (unnumbered) | (roll) | | Biohazard specimen bag (large) | (orde |
| Polymedco Fit Kit | (order as ea) | QTY | MISCELLANEOUS | |
| MICROBIOLOGY | | | Breath test kit | |
| Standard tip Eswab | (order as ea) | | Specimen log book | (orde |
| Nasopharyngeal Eswab | (order as ea) | | Pediatric lead kit | (orde |
| Blood culture bottle AFB | (order as ea) | | Pediatric lead kit (FILTER PAPER ONLY) | (orde |
| Blood culture bottle aerobic (green) | (order as ea) | | 72 hour stool container | |
| Blood culture bottle anaerobic (orange) | (order as ea) | | Fetal fibronectin collection kit | (orde |
| Blood culture prep kit (Chloraprep) | (order as ea) | | Panorama prenatal screen | (orde |
| Blood culture device | (order as ea) | | Horizon collection kit | (orde |
| Blood transfer device | (order as ea) | | Horizon/Panorama combo kit | (orde |
| Hologic aptima swab kit multitest/vaginal #301162 (pink) | (order as ea) | | Cortisol saliva kit | (orde |
| Hologic aptima swab kit unisex #301041 (purple) | (order as ea) | QTY | HISTOLOGY | loide |
| Hologic aptima swab kit unisex #301041 (purple) Hologic aptima urine kit M/F #301040 (yellow) | (order as ea) | VII | Container with formalin 20 mL | (orde |
| | (order as ea) | | Container with formalin 60 mL | |
| Pinworm paddle C & S, O & P stool kit | | | | (orde |
| | (order as ea) | ┨──┤ | Container with formalin 180 mL | (orde |
| STD multi-collection kit (Alinity) | (order as ea) | | Container with formalin 480 mL | (orde |
| UTM viral transport tube | (order as ea) | | Container dry 180 mL | (orde |
| URINE COLLECTION | | | Container dry 32 oz | (orde |
| Sterile orange top urine container - Max Limit 200 | (order as ea) | | Container dry 83 oz | (orde |
| 24 hour urine container - Max Limit 10 | (order as ea) | | Container dry 1 gallon/lids | (orde |
| Urine hat - max limit 25 | (order as ea) | | Flow Media/Hank solution | (orde |
| Urine C & S kit (Gray tube with transfer straw) | (order as ea) | | Michels Solution | (orde |
| Pediatric urine collection kit (PUCK) | (order as ea) | QTY | LABELS | |
| Acetic acid (15 mL) | (6/pkg) | | Specimen label - Epic care link 1" core 300778 | |
| Acetic acid (25 mL) | (6/pkg) | | Manual specimen labels | - |
| Boric acid 4 - 2.5 gram tablets | (order as ea) | | Patient Registry 3x1 - 1" core 300770 | |
| Sodium carbonate 5 grams | (order as ea) | | Specimen label - Epic care link 3" core 300745 | |
| Sodium hydroxide 10 mL | (order as ea) | QTY | OTHER | |
| | (order as ea) | | | |
| Unne vellow vaculainer iu mi | | I | | |
| Urine yellow vacutainer 10 mL | . , | \rightarrow | | |
| Urine transfer straw | (order as ea) | | | |